



红叶 (美国)
RED LEAF USA

美国明尼苏达州圣保罗哈德曼大街南154号, 邮编55075
+1 651 340 1542 Whatsapp +1 952 232 7666
154 Hardman Avenue South Saint Paul,MN 55075,U.S.A

苏州红叶医疗器械有限公司
SUZHOU REDLEAF MEDICAL EQUIPMENT CO.,LTD.

中国江苏省张家港市锦丰镇杨锦公路427号
+86-512-5853 1818
No.427 Yangjin Road,Jinfeng Town,Zhangjiagang City,Jiangsu Province,China.

重庆熙科医疗科技有限公司
CHONGQING XIKE MEDICAL TECHNOLOGY CO.,LTD.

中国重庆市北碚区京东方大道326、324号
+86-23-6826 3062
326 and 324 BOE Road,Beipei District,Chongqing,China.

marketing@chinaredleaf.com
admin@xikemedical.com.cn
www.redleaf.cn



外固定支架 产品手册 EXTERNAL FIXATION PRODUCT CONTENTS

苏州红叶医疗器械有限公司
重庆熙科医疗科技有限公司
SUZHOU REDLEAF MEDICAL EQUIPMENT CO.,LTD.
CHONGQING XIKE MEDICAL TECHNOLOGY CO.,LTD.



公司简介 COMPANY PROFILE



30

成立 30 年
Founded 30 Years

60⁺

国家和地区
Countries & Regions

27,000

覆盖区域
Covered Area

100,000

年生产能力
Annual Production Capacity

苏州红叶医疗器械有限公司成立于 2004 年，是一家由江苏日新医疗设备股份有限公司全资控股的子公司。成立初期为苏州奥萨图有限公司，专业生产病床和医用家具，后转型为专业骨科内植入等高质耗材。像各类耗材厂家一样，经过多年的研发和产品更新迭代，公司取得了欧盟、俄罗斯等国家的准入，例如外固定支架、创伤、脊柱等产品。2019 年，公司取得欧盟 MDD CE 证（包括脊柱后路固定系统、脊柱前路固定系统、椎间融合器、髓内钉、金属股骨颈固定钉、金属骨针、金属接骨螺钉、金属解剖型接骨板、金属空心接骨螺钉、锁定板、不可吸收带线锚钉、外固定支架、颅颌面接骨板钉）。2024 年，又取得了 NMPA 中国医疗器械注册证（灭菌金属骨针、灭菌组合式外固定支架、灭菌单臂一体骨科外固定支架）。

2023 年 11 月，总公司江苏日新医疗设备股份有限公司将重庆熙科医疗器械有限公司全资收购。该公司成立于 2020 年，核心产品围绕足踝病症展开，主要包括有踝关节炎、关节软骨及韧带运动损伤、青少年畸形足等足部矫形、创伤骨折等医疗所需的足踝细分领域的全系列产品，力争成为该领域的旗帜企业。目前，除了常规的数控工制造平台，该公司还开发实现了 3D 打印增材制造，洁净注塑制造，运医编线制造等特殊工艺制造技术，均已实现产品化能力的突破。公司已经陆续取得了包括 12 张三类医疗器械注册证书、2 张二类医疗器械注册证书、2 张一类医疗器械备案在内的 16 张证书。


依托日新医疗的有力支持，红叶医疗和熙科医疗将携手并进，共同建立完整的骨科高端内植入耗材产品体系。

Suzhou Redleaf Medical Equipment Co., Ltd, established in 2004, is a wholly-owned subsidiary of Jiangsu Rixin Medical Equipment Co., Ltd. At the beginning of its establishment, it was named Suzhou Osatu Co., Ltd., specializing in the production of hospital beds and medical furniture, and later transformed into professional orthopedic implants and other high-quality consumables. Like all kinds of consumables manufacturers, after years of research and development and product update iteration, the company has obtained the access of the European Union, Russia and other countries, such as external fixation, trauma, spine and other products. In 2019, the company obtained the EU MDD CE Certificate (including posterior spinal fixation system, anterior spinal fixation system, centrum fusion device, interlocking nail, metal nail of neck fixation of femur, intramedullary pin, metal bone screw, metallic anatomical bone plate, cannulated bone screw, locking plate, non-absorbable suture anchor, orthopedic external fixation system, craniomaxillofacial plate-screw system). In 2024, it also obtained the NMPA China Medical Device Registration Certificate (sterilized metal bone pin, sterilized combined external fixation, sterilized single arm integrated orthopedic external fixation).

In November 2023, Jiangsu Rixin Medical Equipment Co., Ltd. wholly acquired Chongqing Xike Medical Technology Co., Ltd. Founded in 2020, the company's core products focus on foot and ankle diseases, including ankle arthritis, sports injuries of articular cartilage and ligaments, foot orthotics such as teenage deformed feet, trauma fractures and other medical needs of the foot and ankle segments of a full range of products, and strive to become a banner enterprise in this field. At present, in addition to the conventional CNC manufacturing platform, the company has also developed and realized 3D printing additive manufacturing, clean injection molding manufacturing, medical braided wire manufacturing and other special process manufacturing technologies, which have achieved breakthroughs in production. Relying on the strong support of Rixin Medical, Redleaf Medical and Xike Medical will work hand in hand to establish a complete product system of high-end internal implant consumables in orthopedics.

Relying on the strong support of Rixin Medical, Redleaf Medical and Xike Medical will work hand in hand to establish a complete product system of high-end internal implant consumables in orthopedics.





精密设备

PRECISION EQUIPMENT

Arcam EBM
A GE Additive Company

Q10plus



CONTENTS

目录

01	腕关节支架系列 Wrist External Fixator Series
03	微动力支架系列 Dynamic Axial External Fixator Series
04	跨关节支架系列 Transarticular External Fixator Series
05	骨盆支架&骨搬运支架系列 Pelvic & Bone Transport External Fixator Series
06	微型支架系列 Mini External Fixator Series
08	标准支架系列 Standard External Fixator Series
10	组合式外固定支架ZI系列 Combinational External Fixator ZI Series
23	组合式外固定支架ZII系列 Combinational External Fixator ZII Series
27	组合式外固定支架ZIII系列 Combinational External Fixator ZIII Series
32	组合式外固定支架ZIV系列 Combinational External Fixator ZIV Series
35	组合式外固定支架ZV系列 Combinational External Fixator ZV Series
41	组合式外固定支架H系列 Combinational External Fixator H Series
46	钻头 Drill Bits
47	骨牵引针 Bone Screws
50	自攻自钻骨牵引针 Self-tapping self-drilling Bone Screws
52	克氏针 K-wires

腕关节支架系列

WRIST EXTERNAL FIXATOR SERIES

腕关节标准型

Wrist external fixator (Standard type)

产品编号:2201341002160 2201341002161

型号规格:III05 216mm

骨牵引针:Φ2.5, Φ3.5

配套器械: 3.5mm上钉器, 4mm L型扳手

Product No:2201341002160 2201341002161

Model&Specification:III05 216mm

Matching bone screws:Φ2.5, Φ3.5

Matching instruments: 3.5mm screwdriver, 4mm L spanner



腕关节L型

Wrist external fixator (L type)

产品编号:2202341002160 2202341002161

型号规格:III02 216mm

骨牵引针:Φ3.5

配套器械: 3.5mm上钉器, 4mm L型扳手

Product No:2202341002160 2202341002161

Model&Specification:III02 216mm

Matching bone screws:Φ3.5

Matching instruments: 3.5mm screwdriver, 4mm L spanner



腕关节半球型

Wrist external fixator (Hemispherical type)

产品编号:2203341001950 2203341001951

型号规格:III03 195mm

骨牵引针:Φ2.5, Φ3.5


配套器械: 2.5mm上钉器, 3.5mm上钉器, 4mm L型扳手, 10mm开口扳手

Product No:2203341001950 2203341001951

Model&Specification:III03 195mm

Matching bone screws:Φ2.5, Φ3.5

Matching instruments: 2.5mm screwdriver, 3.5mm screwdriver, 4mm L spanner, 10mm Open wrench



腕关节单针型

Wrist external fixator (Single-pin type)

产品编号:2204341002160 2204341002161

型号规格:III06 216mm

骨牵引针:Φ2.5, Φ3.5


配套器械: 3.5mm上钉器, 4mm L型扳手

Product No:2204341002160 2204341002161

Model&Specification:III06 216mm

Matching bone screws:Φ2.5, Φ3.5

Matching instruments: 3.5mm screwdriver, 4mm L spanner



备注: 产品编号末位数0代表非灭菌, 末位数1代表灭菌。

Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

腕关节支架系列

WRIST EXTERNAL FIXATOR SERIES

前臂架

Forearm external fixator

产品编号:2205341001600 2205341001601

型号规格:I02 160mm

骨牵引针:Φ3.5


配套器械: 3.5mm上钉器, 5mm L型扳手, 10mm开口扳手

Product No:2205341001600 2205341001601

Model&Specification:I02 160mm

Matching bone screws:Φ3.5

Matching instruments: 3.5mm screwdriver, 5mm L spanner, 10mm Open wrench



腕关节透光型

Wrist external fixator (Radiolucent type)

产品编号:2206341001500 2206341001501

型号规格:III07 150mm

骨牵引针:Φ2.5, Φ3.5


配套器械: 3.5mm上钉器, 4mm L型扳手

Product No:2206341001500 2206341001501

Model&Specification:III07 150mm

Matching bone screws:Φ2.5, Φ3.5

Matching instruments: 3.5mm screwdriver, 4mm L spanner



腕关节支透光多关节型

Wrist external fixator (Radiolucent multi-joint type)

产品编号:2207341001950 2207341001951

型号规格:III08 195mm

骨牵引针:Φ2.5, Φ3.5


配套器械: 3.5mm上钉器, 5mm L型扳手

Product No:2207341001950 2207341001951

Model&Specification:III08 195mm

Matching bone screws:Φ2.5, Φ3.5

Matching instruments: 3.5mm screwdriver, 5mm L spanner



备注: 产品编号末位数0代表非灭菌, 末位数1代表灭菌。

Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

微动力支架系列
DYNAMIC AXIAL EXTERNAL FIXATOR SERIES

微动支架
Dynamic axial external fixator (Straight type)

产品编号 Product No.	型号规格 Model&Specification
2208341003050 2208341003051	II01 XL 305mm
2208341002650 2208341002651	II01 L 265mm
2208341002420 2208341002421	II01 M 242mm
2208341002080 2208341002081	II01 S 208mm
2208341001750 2208341001751	II01 XS 175mm

骨牵引针:Φ6/Φ5/Φ4
配套器械: 6mm上钉器/5mm上钉器/4mm上钉器; 6mm L型扳手/5mm L型扳手/4mm L型扳手

Matching bone screws:Φ6/Φ5/Φ4
Matching instruments: 6mm screwdriver/5mm screwdriver/4mm screwdriver; 6mm L spanner/5mm L spanner/4mm L spanner



微动支架
Dynamic axial external fixator (T type)

产品编号 Product No.	型号规格 Model&Specification
2209341002680 2209341002681	II01 XL 268mm
2209341002280 2209341002281	II01 L 228mm
2209341001950 2209341001951	II01 M 195mm
2209341001750 2209341001751	II01 S 175mm
2209341001450 2209341001451	II01 XS 145mm

骨牵引针:Φ6/Φ5/Φ4
配套器械: 6mm上钉器/5mm上钉器/4mm上钉器; 6mm L型扳手/5mm L型扳手/4mm L型扳手

Matching bone screws:Φ6/Φ5/Φ4
Matching instruments: 6mm screwdriver/5mm screwdriver/4mm screwdriver; 6mm L spanner/5mm L spanner/4mm L spanner




微动髋关节型
Dynamic axial external fixator (Hip joint type)

产品编号 Product No.	型号规格 Model&Specification
2210341002750 2210341002751	II02 XL 275mm
2210341002350 2210341002351	II02 L 235mm

骨牵引针:Φ6
配套器械: 6mm上钉器, 6mm L型扳手

Matching bone screws:Φ6
Matching instruments: 6mm screwdriver, 6mm L spanner




备注: 产品编号末位数0代表非灭菌, 末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

跨关节支架系列
TRANSARTICULAR EXTERNAL FIXATOR SERIES

粗隆间支架
Petrochanteric external fixator

产品编号:2211341000960 2211341000961
型号规格:II03 96mm
骨牵引针:Φ6
配套器械: 6mm上钉器, 5mm L型扳手
Product No:2211341000960 2211341000961
Model&Specification:II03 96mm
Matching bone screws:Φ6
Matching instruments: 6mm screwdriver, 5mm L spanner



肘关节型
Elbow external fixator

产品编号:2212341003400 2212341003401
型号规格:III10 340mm
骨牵引针:Φ5
配套器械: 5mm上钉器, 5mm L型扳手, 3mm L型扳手
Product No:2212341003400 2212341003401
Model&Specification:III10 340mm
Matching bone screws:Φ5
Matching instruments: 5mm screwdriver, 5mm L spanner, 3mm L spanner



膝关节型
Knee external fixator

产品编号:2213341002700 2213341002701
型号规格:III01 270mm
骨牵引针:Φ6
配套器械: 6mm上钉器, 6mm L型扳手, 3mm L型扳手
Product No:2213341002700 2213341002701
Model&Specification:III01 270mm
Matching bone screws:Φ6
Matching instruments: 6mm screwdriver, 6mm L spanner, 3mm L spanner



踝关节型
Ankle external fixator

产品编号 Product No.	型号规格 Model&Specification
2225341002900 2225341002901	III09 L 290mm
2225341002650 2225341002651	III09 M 265mm

骨牵引针:Φ6/Φ5
配套器械: 6mm上钉器/5mm上钉器; 6mm L型扳手

Matching bone screws:Φ6/Φ5
Matching instruments: 6mm screwdriver/5mm screwdriver; 6mm L spanner



备注: 产品编号末位数0代表非灭菌, 末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

骨盆支架&骨搬运支架系列

PELVIC & BONE TRANSPORT EXTERNAL FIXATOR SERIES

两节型
Pelvic external fixator (Two-section type)

产品编号:2214341002400 2214341002401
型号规格:IV02 240mm
骨牵引针:Φ6
配套器械: 6mm上钉器, 6mm L型扳手

Product No:2214341002400 2214341002401
Model&Specification:IV02 240mm
Matching bone screws:Φ6
Matching instruments: 6mm screwdriver, 6mm L spanner



三节组合型
Pelvic external fixator (Three-section combination type)

产品编号:2215341005200 2215341005201
型号规格:IV03 520mm
骨牵引针:Φ6
配套器械: 6mm上钉器, 6mm L型扳手


Product No:2215341005200 2215341005201
Model&Specification:IV03 520mm
Matching bone screws:Φ6
Matching instruments: 6mm screwdriver, 6mm L spanner



骨搬运支架
Bone transport external fixator

产品编号:2216341003000 2216341003001
型号规格:I04 L 300mm
骨牵引针:Φ6
配套器械: 6mm上钉器, 6mm L型扳手

Product No:2216341003000 2216341003001
Model&Specification:I04 L 300mm
Matching bone screws:Φ6
Matching instruments: 6mm screwdriver, 6mm L spanner



骨搬运支架
Bone transport external fixator

产品编号:2216341002500 2216341002501
型号规格:I04 S 250mm
骨牵引针:Φ4
配套器械: 4mm上钉器, 5mm L型扳手

Product No:2216341002500 2216341002501
Model&Specification:I04 S 250mm
Matching bone screws:Φ4
Matching instruments: 4mm screwdriver, 5mm L spanner



备注: 产品编号末位数0代表非灭菌, 末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

微型支架系列

MINI EXTERNAL FIXATOR SERIES

关节型
Mini external fixator (Articular type)

产品编号:2217341001000 2217341001001
型号规格:V01 100mm
骨牵引针:Φ2
配套器械: 2mm上钉器, 3mm L型扳手

Product No:2217341001000 2217341001001
Model&Specification:V01 100mm
Matching bone screws:Φ2
Matching instruments: 2mm screwdriver, 3mm L spanner



直型
Mini external fixator (Straight type)

产品编号:2218341000800 2218341000801
型号规格:V02 80mm
骨牵引针:Φ2
配套器械: 2mm上钉器, 3mm L型扳手


Product No:2218341000800 2218341000801
Model&Specification:V02 80mm
Matching bone screws:Φ2
Matching instruments: 2mm screwdriver, 3mm L spanner



万向加压延长型
Mini external fixator (Universal compression extension type)

产品编号:2219341000700 2219341000701
型号规格:V03 70mm
骨牵引针:Φ2.5
配套器械: 2mm上钉器, 3mm L型扳手


Product No:2219341000700 2219341000701
Model&Specification:V03 70mm
Matching bone screws:Φ2.5
Matching instruments: 2mm screwdriver, 3mm L spanner



万向型
Mini external fixator (Universal type)

产品编号:2220341000450 2220341000451
型号规格:V04 45mm
骨牵引针:Φ2
配套器械: 2mm上钉器, 3mm L型扳手

Product No:2220341000450 2220341000451
Model&Specification:V04 45mm
Matching bone screws:Φ2
Matching instruments: 2mm screwdriver, 3mm L spanner



备注: 产品编号末位数0代表非灭菌, 末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

微型支架系列
MINI EXTERNAL FIXATOR SERIES

组合型

Mini external fixator (Combined type)

产品编号:2221341000800 2221341000801

型号规格:V05 80mm

骨牵引针:Φ2

配套器械: 2mm上钉器, 3mm L型扳手

Product No:2221341000800 2221341000801

Model&Specification:V05 80mm

Matching bone screws:Φ2

Matching instruments: 2mm screwdriver, 3mm L spanner



轨道直型

Mini external fixator (Track straight type)

产品编号:2222341000720 2222341000721

型号规格:V06 72mm

骨牵引针:Φ2

配套器械: 2mm上钉器, 3mm L型扳手

Product No:2222341000720 2222341000721

Model&Specification:V06 72mm

Matching bone screws:Φ2

Matching instruments: 2mm screwdriver, 3mm L spanner



轨道水平轴型

Mini external fixator (Track horizontal axis type)

产品编号:2223341008850 2223341008851

型号规格:V07 88.5mm

骨牵引针:Φ2

配套器械: 2mm上钉器, 3mm L型扳手

Product No:2223341008850 2223341008851

Model&Specification:V07 88.5mm

Matching bone screws:Φ2

Matching instruments: 2mm screwdriver, 3mm L spanner



轨道垂直轴型

Mini external fixator (Track vertical axis type)

产品编号:2224341006950 2224341006951

型号规格:V08 69.5mm

骨牵引针:Φ2


配套器械: 2mm上钉器, 3mm L型扳手

Product No:2224341006950 2224341006951

Model&Specification:V08 69.5mm

Matching bone screws:Φ2

Matching instruments: 2mm screwdriver, 3mm L spanner



备注: 产品编号末位数0代表非灭菌, 末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

标准支架系列
STANDARD EXTERNAL FIXATOR SERIES

中间夹针

External fixator (Straight pin-in-middle type)

产品编号 Product No.	型号规格 Model&Specification
2226341002850 2226341002851	I06 L 285mm
2226341002800 2226341002801	I06 M 280mm
2226341002520 2226341002521	I06 S 252mm

骨牵引针:Φ6/Φ5

配套器械: 6mm上钉器/5mm上钉器; 6mm L型扳手/5mm L型扳手

Matching bone screws:Φ6/Φ5

Matching instruments: 6mm screwdriver/5mm screwdriver; 6mm L spanner/5mm L spanner



三维型

External fixator (3D type)

产品编号 Product No.	型号规格 Model&Specification
2227341003350 2227341003351	I07 L 335mm
2227341002700 2227341002701	I07 M 270mm
2227341002250 2227341002251	I07 S 225mm
2227341001700 2227341001701	I07 XS 170mm

骨牵引针:Φ6/Φ5/Φ4

配套器械: 6mm上钉器/5mm上钉器/4mm上钉器; 6mm L型扳手/5mm L型扳手/4mm L型扳手; 10-14开口扳手

Matching bone screws:Φ6/Φ5/Φ4

Matching instruments: 6mm screwdriver/5mm screwdriver/4mm screwdriver; 6mm L spanner/5mm L spanner/4mm L spanner; 10-14 open wrench



三节型

External fixator (Three-section straight type)

产品编号 Product No.	型号规格 Model&Specification
2228341003300 2228341003301	I08 L 330mm
2228341003050 2228341003051	I08 M 305mm

骨牵引针:Φ6

配套器械: 6mm上钉器, 6mm L型扳手

Matching bone screws:Φ6

Matching instruments: 6mm screwdriver, 6mm L spanner



三节型

External fixator (Three-section T type)

产品编号 Product No.	型号规格 Model&Specification
2229341002950 2229341002951	II05 L 295mm
2229341002700 2229341002701	II05 M 270mm

骨牵引针:Φ6

配套器械: 6mm上钉器, 6mm L型扳手

Matching bone screws:Φ6

Matching instruments: 6mm screwdriver, 6mm L spanner



备注: 产品编号末位数0代表非灭菌, 末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

标准支架系列

STANDARD EXTERNAL FIXATOR SERIES

关节型

External fixator (Articular type)

产品编号:2230341001900 2230341001901

型号规格:I09 190mm

骨牵引针:Φ6

配套器械: 6mm上钉器, 6mm L型扳手

Product No:2230341001900 2230341001901

Model&Specification:I09 190mm

Matching bone screws:Φ6

Matching instruments: 6mm screwdriver, 6mm L spanner



斜型

External fixator (Oblique type)

产品编号 Product No.	型号规格 Model&Specification
2231341002700 2231341002701	I10 L 270mm
2231341002450 2231341002451	I10 M 245mm



骨牵引针:Φ6

配套器械: 6mm上钉器, 6mm L型扳手

Matching bone screws:Φ6

Matching instruments: 6mm screwdriver, 6mm L spanner

直型

Stainless steel straight external fixator

产品编号 Product No.	型号规格 Model&Specification
223230002500 223230002501	I11 L 250mm
223230002100 223230002101	I11 M 210mm
223230001850 223230001851	I11 S 185mm



骨牵引针:Φ6/Φ5

配套器械: 6mm上钉器/5mm上钉器; 5mm L型扳手/4mm L型扳手; 10-16开口扳手/10-14开口扳手

Matching bone screws:Φ6/Φ5

Matching instruments: 6mm screwdriver/5mm screwdriver; 5mm L spanner/4mm L spanner; 10-16 open wrench/10-14 open wrench

万向型

Stainless steel universal external fixator

产品编号 Product No.	型号规格 Model&Specification
223330002500 223330002501	I12 L 250mm
223330002100 223330002101	I12 M 210mm
223330001850 223330001851	I12 S 185mm



骨牵引针:Φ6/Φ5

配套器械: 6mm上钉器/5mm上钉器; 5mm L型扳手/4mm L型扳手; 10-16开口扳手/10-14开口扳手

Matching bone screws:Φ6/Φ5

Matching instruments: 6mm screwdriver/5mm screwdriver; 5mm L spanner/4mm L spanner; 10-16 open wrench/10-14 open wrench

备注: 产品编号末位0代表非灭菌, 末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZI系列

COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

针杆夹 (单针式)

Pin to rod coupling (Single-pin type)

产品编号 Product No.	型号规格 Model&Specification
2101341005030 2101341005031	ZIP01 Φ5/Φ3
2101341005040 2101341005041	ZIP01 Φ5/Φ4



杆杆夹 (卡入式)

Rod to rod coupling (Snap-in type)

产品编号:2102341005050 2102341005051

型号规格:ZIR01 Φ5/Φ5

Product No:2102341005050 2102341005051

Model&Specification:ZIR01 Φ5/Φ5



针杆夹 (近关节)

Pin to rod coupling (Periarticular type)

产品编号 Product No.	型号规格 Model&Specification
2103341005030 2103341005031	ZIP03 Φ5/Φ3
2103341005040 2103341005041	ZIP03 Φ5/Φ4



针杆夹 (多针式)

Pin to rod coupling (Multi-pin type)

产品编号:2104341000040 2104341000041

型号规格:ZIP02 4holes

Product No:2104341000040 2104341000041

Model&Specification:ZIP02 4holes



连接杆 (支柱)

Connecting rod (Post)

产品编号:210530005000 210530005001

型号规格:ZIL03 Φ5/0°

Product No:210530005000 210530005001

Model&Specification:ZIL03 Φ5/0°



备注: 产品编号末位0代表非灭菌, 末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

连接杆（支柱）
Connecting rod (Post)

产品编号:210530005300 210530005301
型号规格:ZIL03 Φ5/30°

Product No:210530005300 210530005301
Model&Specification:ZIL03 Φ5/30°



连接杆（直杆）
Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
210641005100 ~ 2106410052000	ZIL02 Φ5/10 ~ 200mm
210641005101 ~ 2106410052001	(at intervals of 5mm)



连接杆（直杆）
Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
210670005100 ~ 2106700052000	ZIL02 Φ5/10 ~ 200mm
210670005101 ~ 2106700052001	(at intervals of 5mm)



跨关节连接件（肘关节）
Transarticular connector (Elbow)

产品编号:2107300000050 2107300000051
型号规格:ZVIIG03 Φ5

Product No:2107300000050 2107300000051
Model&Specification:ZVIIG03 Φ5



跨关节连接件（腕关节）
Transarticular connector (Wrist)

产品编号:2108300000050 2108300000051
型号规格:ZVIIG04 Φ5

Product No:2108300000050 2108300000051
Model&Specification:ZVIIG04 Φ5



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

杆杆夹（卡入式）
Rod to rod coupling (Snap-in type)

产品编号:2109341015050 2109341015051
型号规格:ZIR01 Φ15/Φ5

Product No:2109341015050 2109341015051
Model&Specification:ZIR01 Φ15/Φ5



连接杆（加压延长杆单杆）
Connecting rod (Compression extension single rod)

产品编号:2110341000150 2110341000151
型号规格:ZVL05 Φ15

Product No:2110341000150 2110341000151
Model&Specification:ZVL05 Φ15



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES


尺桡骨跨关节支架
Ulna & radius transarticular external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ5/Φ5	4
Pin to rod coupling (Multi-pin type) ZIP02 4holes	2
Connecting rod (Post) ZIL03 Φ5/30°	4
Connecting rod (Straight rod) ZIL02 Φ5	2
Bone screw Φ3	4




前臂支架
Forearm external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ5/Φ5	2
Pin to rod coupling (Multi-pin type) ZIP02 4holes	2
Connecting rod (Post) ZIL03 Φ5/30°	2
Connecting rod (Straight rod) ZIL02 Φ5	1
Bone screw Φ3	4




尺桡骨支架
Ulna & radius external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ5/Φ5	1
Pin to rod coupling (Single-pin type) ZIP01 Φ5/Φ3	2
Pin to rod coupling (Periarticular type) ZIP03 Φ5/Φ3	1
Connecting rod (Straight rod) ZIL02 Φ5	1
Bone screw Φ3	4



腕关节支架
Wrist joint external fixator

配置 Configuration	数量 Qty
Transarticular connector (Wrist) ZVIIIG04 Φ5	1
Pin to rod coupling (Single-pin type) ZIP01 Φ5/Φ3	4
Self-tapping self-drilling bone screw Φ3	4



组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

肘关节支架
Elbow joint external fixator

配置 Configuration	数量 Qty
Transarticular connector (Elbow) ZVIIIG03 Φ5	1
Pin to rod coupling (Single-pin type) ZIP01 Φ5/Φ3	4
Self-tapping self-drilling bone screw Φ4	4



肱骨干支架
Humeral shaft external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ5/Φ5	2
Pin to rod coupling (Multi-pin type) ZIP02 4holes	2
Connecting rod (Post) ZIL03 Φ5/30°	2
Connecting rod (Straight rod) ZIL02 Φ5	1
Bone screw Φ4	4




肱骨近端支架
Proximal humeral external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ5/Φ5	4
Pin to rod coupling (Single-pin type) ZIP01 Φ5/Φ3	2
Pin to rod coupling (Multi-pin type) ZIP02 4holes	1
Connecting rod (Straight rod) ZIL02 Φ5	3
Connecting rod (Post) ZIL03 Φ5/30°	2
Self-tapping self-drilling bone screw Φ4	4




肱骨支架
Humeral external fixator

配置 Configuration	数量 Qty
Connecting rod (Compression extension rod) ZVL05 Φ15	1
Rod to rod coupling (Snap-in type) ZIR01 Φ15/Φ5	2
Pin to rod coupling (Multi-pin type) ZIP02 4holes	2
Connecting rod (Post) ZIL03 Φ5/30°	2
Bone screw Φ4	4



组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

组合式外固定ZI系列器械包I（标配）
QE100 Combinational External Fixation ZI Instrument Kit I（Standard configuration）




序号 No.	产品编号 Product No.	产品名称 Product Name		型号规格 Specification	数量 Qty
1	2102341005050	杆杆夹（卡入式）	Rod to rod coupling (Snap-in type)	ZIR01 Φ5/Φ5	12
2	2101341005030	针杆夹（单针式）	Pin to rod coupling (Single-pin type)	ZIP01 Φ5/Φ3	12
3	2104341000040	针杆夹（多针式）	Pin to rod coupling (Multi-pin type)	ZIP02 4holes	3
4	2103341005030	针杆夹（近关节）	Pin to rod coupling (Periarticular type)	ZIP03 Φ5/Φ3	1
5	210530005000	连接杆（支柱）	Connecting rod (Post)	ZIL03 Φ5/0°	2
6	210530005300	连接杆（支柱）	Connecting rod (Post)	ZIL03 Φ5/30°	4
7	G23080S	自攻自钻骨牵引针	Self-tapping self-drilling bone screw	Φ3.0×80	4
	G230100S	自攻自钻骨牵引针	Self-tapping self-drilling bone screw	Φ3.0×100	4
	G24080S	自攻自钻骨牵引针	Self-tapping self-drilling bone screw	Φ4.0×80	4
	G240100S	自攻自钻骨牵引针	Self-tapping self-drilling bone screw	Φ4.0×100	4
	G240120S	自攻自钻骨牵引针	Self-tapping self-drilling bone screw	Φ4.0×120	4
8	QE100-01	定位器	Guider	Φ3	1
	QE100-02	定位器	Guider	Φ4	1
9	2106700051200	连接杆（直杆）	Connecting rod (Straight rod)	ZIL02 Φ5/120mm	4
	2106700051500	连接杆（直杆）	Connecting rod (Straight rod)	ZIL02 Φ5/150mm	4
	2106700051800	连接杆（直杆）	Connecting rod (Straight rod)	ZIL02 Φ5/180mm	2
	2106700052000	连接杆（直杆）	Connecting rod (Straight rod)	ZIL02 Φ5/200mm	4
10	QE100-03	T型扳手	T wrench	#5	1
11	QE100-04	稳定复位扳手	Stablization reduction wrench	#15	1
12	210730000050	跨关节连接件（肘关节）	Transarticular connector (Elbow)	ZVIG03 Φ5	1
13	QE100-05	骨圆针连接头	Adaptor	Φ3	1
	QE100-06	骨圆针连接头	Adaptor	Φ4	1
14	QE100-07	上钉器	Screwdriver	Φ3	1
	QE100-08	上钉器	Screwdriver	Φ4	1
15	QE100-09	弓形钻架	Bow-shaped frame	/	1
16	QE100-10	拇指轮	Thumb wheel	#5/7	1
17	QE100-11	开口扳手	Open spanner	#5/7	1
18	QE100-12	器械盒	Instrument box	/	1

组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

针杆夹（单针式）
Pin to rod coupling (Single-pin type)

产品编号 Product No.	型号规格 Model&Specification
2101341003150 2101341003151	ZIP01 Φ3/Φ1.5
2101341003020 2101341003021	ZIP01 Φ3/Φ2



连接杆（直杆）
Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
210641003100 ~ 2106410031000	ZIL02 Φ3/10 ~ 100mm
210641003101 ~ 2106410031001	(at intervals of 5mm)



连接杆（直杆）
Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
210670003100 ~ 2106700031000	ZIL02 Φ3/10 ~ 100mm
210670003101 ~ 2106700031001	(at intervals of 5mm)



指掌骨支架
Mini external fixator

配置 Configuration	数量 Qty
Pin to rod coupling (Single-pin type) ZIP01 Φ3/Φ2	2
Connecting rod (Straight rod) ZIL02 Φ3	1
Bone screw Φ2	4



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

针杆夹（单针式）
Pin to rod coupling (Single-pin type)

产品编号 Product No.	型号规格 Model&Specification
2101341008040 2101341008041	ZIP01 Φ8/Φ4
2101341008050 2101341008051	ZIP01 Φ8/Φ5



杆杆夹（卡入式）
Rod to rod coupling (Snap-in type)

产品编号:2102341008080 2102341008081
型号规格:ZIR01 Φ8/Φ8

Product No:2102341008080 2102341008081
Model&Specification:ZIR01 Φ8/Φ8



针杆夹（多针式）
Pin to rod coupling (Multi-pin type)

产品编号:2104341000050 2104341000051
型号规格:ZIP02 5holes
Product No:2104341000050 2104341000051
Model&Specification:ZIP02 5holes



连接杆（支柱）
Connecting rod (Post)

产品编号:210530008000 210530008001
型号规格:ZIL03 Φ8/0°
Product No:210530008000 210530008001
Model&Specification:ZIL03 Φ8/0°



连接杆（支柱）
Connecting rod (Post)

产品编号:210530008300 210530008301
型号规格:ZIL03 Φ8/30°
Product No:210530008300 210530008301
Model&Specification:ZIL03 Φ8/30°



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

连接杆（直杆）
Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
210641008100 ~ 2106410084000	ZIL02 Φ8/10 ~ 400mm
210641008101 ~ 2106410084001	(at intervals of 5mm)



连接杆（直杆）
Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
210670008100 ~ 2106700084000	ZIL02 Φ8/10 ~ 400mm
210670008101 ~ 2106700084001	(at intervals of 5mm)



跨关节连接件（踝关节）
Transarticular connector (Ankle)

产品编号:211130000080 211130000081
型号规格:ZVIIG02 Φ8
Product No:211130000080 211130000081
Model&Specification:ZVIIG02 Φ8



跨关节连接件（膝关节）
Transarticular connector (Knee)

产品编号 Product No.	型号规格 Model&Specification
211230100080 211230100081	ZVIIG01 Φ8 L
211230200080 211230200081	ZVIIG01 Φ8 R



杆杆夹（卡入式）
Rod to rod coupling (Snap-in type)

产品编号:2109341020080 2109341020081
型号规格:ZIR01 Φ20/Φ8
Product No:2109341020080 2109341020081
Model&Specification:ZIR01 Φ20/Φ8



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

连接杆（加压延长杆单杆）
Connecting rod (Compression extension single rod)

产品编号:2110341000200 2110341000201
型号规格:ZVL05 Φ20
Product No:2110341000200 2110341000201
Model&Specification:ZVL05 Φ20



连接杆（加压延长杆单杆）
Connecting rod (Compression extension single rod)

产品编号:2113341000080 2113341000081
型号规格:ZVL05 Φ8
Product No:2113341000080 2113341000081
Model&Specification:ZVL05 Φ8



连接杆（Φ8半环连接杆）
Connecting rod (Φ8 semi-ring)

产品编号 Product No.	型号规格 Model&Specification
211441081600 211441081601	ZIL01 Φ8×160mm
211441081800 211441081801	ZIL01 Φ8×180mm
211441082000 211441082001	ZIL01 Φ8×200mm



胫骨近端支架
Proximal tibia external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ8/Φ8	4
Pin to rod coupling (Single-pin type) ZIP01 Φ8/Φ5	3
Pin to rod coupling (Multi-pin type) ZIP02 5holes	1
Connecting rod (Post) ZIL03 Φ8/30°	2
Connecting rod (Φ8 semi-ring)	1
Connecting rod (Straight rod) ZIL02 Φ8	2
Self-tapping self-drilling bone screw Φ5	5



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

股骨干支架
Femoral shaft external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ8/Φ8	2
Pin to rod coupling (Single-pin type) ZIP01 Φ8/Φ5	4
Connecting rod (Straight rod) ZIL02 Φ8	3
Self-tapping self-drilling bone screw Φ5	4



标准半环型支架
Standard semi-ring external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ8/Φ8	6
Pin to rod coupling (Multi-pin type) ZIP02 5holes	2
Connecting rod (Post) ZIL03 Φ8/30°	4
Connecting rod (Straight rod) ZIL02 Φ8	3
Self-tapping self-drilling bone screw Φ5	4



骨盆支架
Pelvic external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ8/Φ8	6
Pin to rod coupling (Multi-pin type) ZIP02 5holes	2
Connecting rod (Post) ZIL03 Φ8/30°	2
Connecting rod (Straight rod) ZIL02 Φ8	4
Bone screw Φ5	4



踝关节支架
Ankle joint external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ8/Φ8	3
Pin to rod coupling (Single-pin type) ZIP01 Φ8/Φ5	4
Connecting rod (Straight rod) ZIL02 Φ8	3
Self-tapping self-drilling bone screw Φ5	4



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

膝关节支架
Knee joint external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ8/Φ8	2
Pin to rod coupling (Single-pin type) ZIP01 Φ8/Φ5	4
Transarticular connector (Knee) ZVIG01 Φ8 L	1
Pin to rod coupling (Multi-pin type) ZIP02 5holes	1
Connecting rod (Post) ZIL03 Φ8/30°	1
Connecting rod (Straight rod) ZIL02 Φ8	1
Self-tapping self-drilling bone screw Φ5	6



踝关节融合支架
Ankle fusion external fixator

配置 Configuration	数量 Qty
Pin to rod coupling (Single-pin type) ZIP01 Φ5/Φ3	1
Pin to rod coupling (Single-pin type) ZIP01 Φ8/Φ5	4
Transarticular connector (Ankle) ZVIG02 Φ8	1
Self-tapping self-drilling bone screw Φ5	4
Self-tapping self-drilling bone screw Φ3	1



微动加压支架
Micro-motion compression external fixator


配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIR01 Φ20/Φ8	2
Connecting rod (Compression extension rod) ZVL05 Φ20	1
Pin to rod coupling (Multi-pin type) ZIP02 5holes	2
Connecting rod (Post) ZIL03 Φ8/30°	2
Bone screw Φ5	4



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZI系列
COMBINATIONAL EXTERNAL FIXATOR ZI SERIES

组合式外固定ZI系列器械包II（标配）
QE200 Combinational External Fixation ZI Instrument Kit II（Standard configuration）



序号 No.	产品编号 Product No.	产品名称 Product Name		型号规格 Specification	数量 Qty
1	2102341008080	杆杆夹（卡入式）	Rod to rod coupling (Snap-in type)	ZIR01 Φ8/Φ8	12
2	2101341008050	针杆夹（单针式）	Pin to rod coupling (Single-pin type)	ZIP01 Φ8/Φ5	12
3	2104341000050	针杆夹（多针式）	Pin to rod coupling (Multi-pin type)	ZIP02 5holes	3
4	210530008000	连接杆（支柱）	Connecting rod (Post)	ZIL03 Φ8/0°	2
5	210530008300	连接杆（支柱）	Connecting rod (Post)	ZIL03 Φ8/30°	4
6	G250120S	自攻自钻骨牵引针	Self-tapping self-drilling bone screw	Φ5.0×120	4
	G250150S	自攻自钻骨牵引针	Self-tapping self-drilling bone screw	Φ5.0×150	8
	G250180S	自攻自钻骨牵引针	Self-tapping self-drilling bone screw	Φ5.0×180	4
7	QE200-01	定位器	Guider	Φ5	2
8	2106700082000	连接杆（直杆）	Connecting rod (Straight rod)	ZIL02 Φ8/200mm	4
	2106700082500	连接杆（直杆）	Connecting rod (Straight rod)	ZIL02 Φ8/250mm	4
	2106700083000	连接杆（直杆）	Connecting rod (Straight rod)	ZIL02 Φ8/300mm	4
9	211441081800	连接杆（Φ8半环连接杆）	Connecting rod (Φ8 semi-ring)	ZIL01 Φ8×180mm	2
	211441082000	连接杆（Φ8半环连接杆）	Connecting rod (Φ8 semi-ring)	ZIL01 Φ8×200mm	2
10	QE200-02	稳定复位扳手	Stablization reduction wrench	#22	1
11	QE200-03	扭力扳手	Torque wrench	#7	1
12	QE200-04	拇指轮	Thumb wheel	#5/7	1
13	QE200-05	骨圆针连接头	Adaptor	Φ5	1
14	QE200-06	弓形钻架	Bow-shaped frame	/	1
15	QE200-07	上钉器	Screwdriver	Φ5	1
16	QE200-08	T型扳手	T wrench	#7	1
17	QE200-09	开口扳手	Open spanner	#5/7	1
18	QE200-10	器械盒	Instrument box	/	1

组合式外固定支架ZII系列
COMBINATIONAL EXTERNAL FIXATOR ZII SERIES

针杆夹（单针式）
Pin to rod coupling (Single-pin type)

产品编号 Product No.	型号规格 Model&Specification
2115341008030 2115341008031	ZIIP01 Φ8/Φ3
2115341008040 2115341008041	ZIIP01 Φ8/Φ4
2115341011050 2115341011051	ZIIP01 Φ11/Φ5
2115341011060 2115341011061	ZIIP01 Φ11/Φ6



杆杆夹（穿针式）
Rod to rod coupling (Threading type)

产品编号 Product No.	型号规格 Model&Specification
2116341008080 2116341008081	ZIIR01 Φ8/Φ8
2116341011110 2116341011111	ZIIR01 Φ11/Φ11



针杆夹（双针式）
Pin to rod coupling (Double-pin type)

产品编号 Product No.	型号规格 Model&Specification
2117341008040 2117341008041	ZIIP02 Φ8/Φ4
2117341011050 2117341011051	ZIIP02 Φ11/Φ5



连接杆（直杆）
Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
211841008100 ~ 2118410083000	ZIIL01 Φ8/10 ~ 300mm (at intervals of 5mm)
211841011100 ~ 2118410114000	ZIIL01 Φ11/10 ~ 400mm (at intervals of 5mm)



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZII系列
COMBINATIONAL EXTERNAL FIXATOR ZII SERIES

连接杆（直杆）
Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
211870008100 ~ 2118700083000	ZIIL01 Φ8/10 ~ 300mm (at intervals of 5mm)
211870011100 ~ 2118700114000	ZIIL01 Φ11/10 ~ 400mm (at intervals of 5mm)



连接杆（Φ8半环连接杆）
Connecting rod (Φ8 semi-ring)

产品编号 Product No.	型号规格 Model&Specification
211941081600 211941081601	ZIIL02 Φ8×160mm
211941081800 211941081801	ZIIL02 Φ8×180mm
211941082000 211941082001	ZIIL02 Φ8×200mm



连接杆（弧型杆）
Connecting rod (Arc rod)

产品编号:212041082850 212041082851
型号规格:ZIIL03 Φ8×285mm

Product No:212041082850 212041082851
Model&Specification:ZIIL03 Φ8×285mm



连接杆（梯形杆）
Connecting rod (Trapezoid rod)

产品编号:212141081800 212141081801
型号规格:ZIIL04 Φ8×180mm

Product No:212141081800 212141081801
Model&Specification:ZIIL04 Φ8×180mm



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZII系列
COMBINATIONAL EXTERNAL FIXATOR ZII SERIES

T型支架
T-shaped external fixator

配置 Configuration	数量 Qty
Pin to rod coupling (Single-pin type) ZIIP01 Φ11/Φ5	4
Pin to rod coupling (Double-pin type) ZIIP02 Φ11/Φ5	1
Connecting rod (Straight rod) ZIIL01 Φ11	1
Bone screw Φ5	6



Z型支架
Z-shaped external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Threading type) ZIIR01 Φ11/Φ11	2
Pin to rod coupling (Single-pin type) ZIIP01 Φ11/Φ5	4
Connecting rod (Straight rod) ZIIL01 Φ11	3
Bone screw Φ5	4



直型支架
Straight external fixator


配置 Configuration	数量 Qty
Pin to rod coupling (Single-pin type) ZIIP01 Φ11/Φ5	4
Connecting rod (Straight rod) ZIIL01 Φ11	1
Bone screw Φ5	4



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZII系列
COMBINATIONAL EXTERNAL FIXATOR ZII SERIES

组合式外固定ZII系列器械包（标配）
QE300 Combinational External Fixation ZII Instrument Kit (Standard configuration)




序号 No.	产品编号 Product No.	产品名称 Product Name	型号规格 Specification	数量 Qty
1	2115341008030	针杆夹（单针式）	Pin to rod coupling (Single-pin type)	8
2	2116341008080	杆杆夹（穿针式）	Rod to rod coupling (Threading type)	4
3	2116341011110	杆杆夹（穿针式）	Rod to rod coupling (Threading type)	4
4	2115341011050	针杆夹（单针式）	Pin to rod coupling (Single-pin type)	12
5	2117341011050	针杆夹（双针式）	Pin to rod coupling (Double-pin type)	2
6	G14090S	骨牵引针	Bone screw	4
	G150110S	骨牵引针	Bone screw	4
	G150130S	骨牵引针	Bone screw	12
	G160150S	骨牵引针	Bone screw	4
7	Z28200	骨钻	Drill bit	2
	Z38200	骨钻	Drill bit	2
	Z48200	骨钻	Drill bit	2
8	QE300-01	上钉器	Screwdriver	1
	QE300-02	上钉器	Screwdriver	1
9	QE300-03	导向保护套	Guide sleeve	1
	QE300-04	导向保护套	Guide sleeve	1
	QE300-05	导向保护套	Guide sleeve	1
10	QE300-07	开口扳手	Open spanner	2
11	QE300-06	加压器	Pressurizer	1
12	2118700081500	连接杆（直杆）	Connecting rod (Straight rod)	2
	2118700082000	连接杆（直杆）	Connecting rod (Straight rod)	2
	2118700082500	连接杆（直杆）	Connecting rod (Straight rod)	2
	2118700083000	连接杆（直杆）	Connecting rod (Straight rod)	2
	2118700111500	连接杆（直杆）	Connecting rod (Straight rod)	2
	2118700112000	连接杆（直杆）	Connecting rod (Straight rod)	2
	2118700112500	连接杆（直杆）	Connecting rod (Straight rod)	2
	2118700113000	连接杆（直杆）	Connecting rod (Straight rod)	2
13	2118700113500	连接杆（直杆）	Connecting rod (Straight rod)	2
	QE300-08	器械盒	Instrument box	1

组合式外固定支架ZIII系列
COMBINATIONAL EXTERNAL FIXATOR ZIII SERIES


针杆夹（开口）
Pin to rod coupling (Ringent type)

产品编号 Product No.	型号规格 Model&Specification
2122341008030 2122341008031	ZIIIP01 Φ8/Φ3
2122341008040 2122341008041	ZIIIP01 Φ8/Φ4



针杆夹（开口）
Pin to rod coupling (Ringent type)

产品编号 Product No.	型号规格 Model&Specification
2122341011050 2122341011051	ZIIIP01 Φ11/Φ5
2122341011060 2122341011061	ZIIIP01 Φ11/Φ6



杆杆夹（开口）
Rod to rod coupling (Ringent type)

产品编号 Product No.	型号规格 Model&Specification
2123341008080 2123341008081	ZIIIR01 Φ8/Φ8
2123341011110 2123341011111	ZIIIR01 Φ11/Φ11



针杆夹（双针式）
Pin to rod coupling (Double-pin type)

产品编号 Product No.	型号规格 Model&Specification
2117341008040 2117341008041	ZIIP02 Φ8/Φ4
2117341011050 2117341011051	ZIIP02 Φ11/Φ5



连接杆（直杆）
Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
212441008100 ~ 2124410083000	ZIIIL01 Φ8/10 ~ 300mm (at intervals of 5mm)
212441011100 ~ 2124410114000	ZIIIL01 Φ11/10 ~ 400mm (at intervals of 5mm)



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZIII系列
COMBINATIONAL EXTERNAL FIXATOR ZIII SERIES

连接杆（直杆）
Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
212470008100 ~ 2124700083000	ZIIIL01 Φ8/10 ~ 300mm (at intervals of 5mm)
212470011100 ~ 2124700114000	ZIIIL01 Φ11/10 ~ 400mm (at intervals of 5mm)



连接杆（Φ8半环连接杆）
Connecting rod (Φ8 semi-ring)

产品编号 Product No.	型号规格 Model&Specification
212541081600 212541081601	ZIIIL02 Φ8×160mm
212541081800 212541081801	ZIIIL02 Φ8×180mm
212541082000 212541082001	ZIIIL02 Φ8×200mm



连接杆（弧型杆）
Connecting rod (Arc rod)

产品编号:212641082850 212641082851
型号规格:ZIIIL03 Φ8×285mm

Product No:212641082850 212641082851
Model&Specification:ZIIIL03 Φ8×285mm



连接杆（梯形杆）
Connecting rod (Trapezoid rod)

产品编号:212741081800 212741081801
型号规格:ZIIIL04 Φ8×180mm

Product No:212741081800 212741081801
Model&Specification:ZIIIL04 Φ8×180mm



连接管（加压延长管）
Connecting tube (Compression extension tube)

产品编号:212830083300 212830083301
型号规格:ZVL04 Φ8×330mm

Product No:212830083300 212830083301
Model&Specification:ZVL04 Φ8×330mm



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZIII系列
COMBINATIONAL EXTERNAL FIXATOR ZIII SERIES

半环型支架
Semi-ring external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Ringent type) ZIIIR01 Φ8/Φ8	6
Pin to rod coupling (Ringent type) ZIIIP01 Φ8/Φ4	6
Connecting rod (Φ8 semi-ring) ZIIIL02 Φ8	2
Connecting rod (Straight rod) ZIIIL01 Φ8	1
Connecting tube (Compression extension tube) ZVL04 Φ8	2
Bone screw Φ4	6



髌上支架
Condylar external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Ringent type) ZIIIR01 Φ8/Φ8	3
Pin to rod coupling (Ringent type) ZIIIP01 Φ8/Φ4	6
Connecting rod (Φ8 semi-ring) ZIIIL02 Φ8	1
Connecting rod (Trapezoid rod) ZIIIL04 Φ8	1
Connecting rod (Straight rod) ZIIIL01 Φ8	1
Bone screw Φ4	6



胫骨关节支架
Tibia joint external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Ringent type) ZIIIR01 Φ8/Φ8	3
Pin to rod coupling (Ringent type) ZIIIP01 Φ8/Φ4	7
Connecting rod (Φ8 semi-ring) ZIIIL02 Φ8	1
Connecting rod (Arc rod) ZIIIL03 Φ8	1
Connecting rod (Straight rod) ZIIIL01 Φ8	1
Bone screw Φ4	7



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZIII系列
COMBINATIONAL EXTERNAL FIXATOR ZIII SERIES

Z型（T型）支架
Z-shaped external fixator (T type)

配置 Configuration	数量 Qty
Rod to rod coupling (Ringent type) ZIIIR01 Φ8/Φ8	2
Pin to rod coupling (Ringent type) ZIIIP01 Φ8/Φ4	4
Pin to rod coupling (Double-pin type) ZIIP02 Φ8/Φ4	1
Connecting rod (Straight rod) ZIIIL01 Φ8	3
Bone screw Φ4	6



骨盆支架
Pelvic external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Ringent type) ZIIIR01 Φ11/Φ11	2
Pin to rod coupling (Ringent type) ZIIIP01 Φ11/Φ5	4
Pin to rod coupling (Ringent type) ZIIIP01 Φ8/Φ4	2
Connecting rod (Straight rod) ZIIIL01 Φ11	3
Connecting tube (Compression extension tube) ZVL04 Φ8	1
Bone screw Φ5	4




备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定ZIII系列器械包（标配）

COMBINATIONAL EXTERNAL FIXATOR ZIII SERIES

组合式外固定ZIII系列器械包（标配）

QE400 Combinational External Fixation ZIII Instrument Kit (Standard configuration)



序号 No.	产品编号 Product No.	产品名称 Product Name	型号规格 Specification	数量 Qty
1	2123341008080	杆杆夹（开口）	Rod to rod coupling (Ringent type)	4
2	2122341008040	针杆夹（开口）	Pin to rod coupling (Ringent type)	8
3	2123341011110	杆杆夹（开口）	Rod to rod coupling (Ringent type)	4
4	2122341011050	针杆夹（开口）	Pin to rod coupling (Ringent type)	12
5	2117341011050	针杆夹（双针式）	Pin to rod coupling (Double-pin type)	2
6	G14090S	骨牵引针	Bone screw	4
	G150110S	骨牵引针	Bone screw	4
	G150130S	骨牵引针	Bone screw	6
	G150150S	骨牵引针	Bone screw	2
	G160130S	骨牵引针	Bone screw	6
	G160150S	骨牵引针	Bone screw	2
7	Z28200	骨钻	Drill bit	2
	Z38200	骨钻	Drill bit	2
	Z48200	骨钻	Drill bit	2
8	QE400-01	上钉器	Screwdriver	1
	QE400-02	上钉器	Screwdriver	1
9	QE400-03	导向保护套	Guide sleeve	1
	QE400-04	导向保护套	Guide sleeve	1
	QE400-05	导向保护套	Guide sleeve	1
10	QE400-06	开口扳手	Open spanner	2
11	QE400-07	加压器	Pressurizer	1
12	2124700081500	连接杆（直杆）	Connecting rod (Straight rod)	2
	2124700082000	连接杆（直杆）	Connecting rod (Straight rod)	2
	2124700082500	连接杆（直杆）	Connecting rod (Straight rod)	2
	2124700083000	连接杆（直杆）	Connecting rod (Straight rod)	2
	2124700111500	连接杆（直杆）	Connecting rod (Straight rod)	2
	2124700112000	连接杆（直杆）	Connecting rod (Straight rod)	2
	2124700112500	连接杆（直杆）	Connecting rod (Straight rod)	2
	2124700113000	连接杆（直杆）	Connecting rod (Straight rod)	2
	2124700113500	连接杆（直杆）	Connecting rod (Straight rod)	2
	2124700113500	连接杆（直杆）	Connecting rod (Straight rod)	2
13	QE400-08	器械盒	Instrument box	1

组合式外固定支架ZIV系列

COMBINATIONAL EXTERNAL FIXATOR ZIV SERIES

针杆夹（卡入式）

Pin to rod coupling (Snap-in type)

产品编号 Product No.	型号规格 Model&Specification
2129341005150 2129341005151	ZIVP01 Φ5/Φ1.5
2129341005020 2129341005021	ZIVP01 Φ5/Φ2
2129341008030 2129341008031	ZIVP01 Φ8/Φ3
2129341008040 2129341008041	ZIVP01 Φ8/Φ4
2129341011050 2129341011051	ZIVP01 Φ11/Φ5
2129341011060 2129341011061	ZIVP01 Φ11/Φ6



杆杆夹（卡入式）

Rod to rod coupling (Snap-in type)

产品编号 Product No.	型号规格 Model&Specification
2130341005050 2130341005051	ZIVR01 Φ5/Φ5
2130341008080 2130341008081	ZIVR01 Φ8/Φ8
2130341011110 2130341011111	ZIVR01 Φ11/Φ11



针杆夹（双针式）

Pin to rod coupling (Double-pin type)

产品编号 Product No.	型号规格 Model&Specification
2117341008040 2117341008041	ZIIP02 Φ8/Φ4
2117341011050 2117341011051	ZIIP02 Φ11/Φ5



连接杆（直杆）

Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
213141008100 ~ 2131410083000	ZIVL01 Φ8/10 ~ 300mm
213141008101 ~ 2131410083001	(at intervals of 5mm)
213141011100 ~ 2131410114000	ZIVL01 Φ11/10 ~ 400mm
213141011101 ~ 2131410114001	(at intervals of 5mm)



连接杆（直杆）

Connecting rod (Straight rod)

产品编号 Product No.	型号规格 Model&Specification
213170008100 ~ 2131700083000	ZIVL01 Φ8/10 ~ 300mm
213170008101 ~ 2131700083001	(at intervals of 5mm)
213170011100 ~ 2131700114000	ZIVL01 Φ11/10 ~ 400mm
213170011101 ~ 2131700114001	(at intervals of 5mm)



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZIV系列
COMBINATIONAL EXTERNAL FIXATOR ZIV SERIES

连接杆（弧型杆）
Connecting rod (Arc rod)

产品编号:213341082850 213341082851
型号规格:ZIVL03 Φ8×285mm

Product No:213341082850 213341082851
Model&Specification:ZIVL03 Φ8×285mm



连接杆（梯形杆）
Connecting rod (Trapezoid rod)

产品编号:213441081800 213441081801
型号规格:ZIVL04 Φ8×180mm

Product No:213441081800 213441081801
Model&Specification:ZIVL04 Φ8×180mm



T型支架
T-shaped external fixator

配置 Configuration	数量 Qty
Pin to rod coupling (Snap-in type) ZIVP01 Φ11/Φ5	4
Pin to rod coupling (Double-pin type) ZIIP02 Φ11/Φ5	1
Connecting rod (Straight rod) ZIVL01 Φ11	1
Bone screw Φ5	6



直型架
Straight external fixator

配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIVR01 Φ11/Φ11	1
Pin to rod coupling (Double-pin type) ZIIP02 Φ11/Φ5	4
Connecting rod (Straight rod) ZIVL01 Φ11	2
Bone screw Φ5	4



Z型支架
Z-shaped external fixator


配置 Configuration	数量 Qty
Rod to rod coupling (Snap-in type) ZIVR01 Φ11/Φ11	2
Pin to rod coupling (Double-pin type) ZIIP02 Φ11/Φ5	4
Connecting rod (Straight rod) ZIVL01 Φ11	3
Bone screw Φ5	4



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZIV系列
COMBINATIONAL EXTERNAL FIXATOR ZIV SERIES

组合式外固定ZIV系列器械包（标配）
QE500 Combinational External Fixation ZIV Instrument Kit (Standard configuration)




序号 No.	产品编号 Product No.	产品名称 Product Name		型号规格 Specification	数量 Qty
1	2129341011050	针杆夹（卡入式）	Pin to rod coupling (Snap-in type)	ZIVP01 Φ11/Φ5	12
2	2130341011110	杆杆夹（卡入式）	Rod to rod coupling (Snap-in type)	ZIVR01 Φ11/Φ11	9
3	2117341011050	针杆夹（双针式）	Pin to rod coupling (Double-pin type)	ZIIP02 Φ11/Φ5	2
4	G150110S	骨牵引针	Bone screw	Φ5.0×110	4
	G150130S	骨牵引针	Bone screw	Φ5.0×130	6
	G150150S	骨牵引针	Bone screw	Φ5.0×150	2
	G160130S	骨牵引针	Bone screw	Φ6.0×130	6
	G160150S	骨牵引针	Bone screw	Φ6.0×150	2
5	Z38200	骨钻	Drill bit	Φ3.8×200	2
	Z48200	骨钻	Drill bit	Φ4.8×200	2
6	QE500-01	上钉器	Screwdriver	Φ5	1
	QE500-02	上钉器	Screwdriver	Φ6	1
7	QE500-03	导向保护套	Guide sleeve	Φ4	1
	QE500-04	导向保护套	Guide sleeve	Φ5	1
	QE500-05	导向保护套	Guide sleeve	Φ6	1
8	QE500-06	开口扳手	Open spanner	#10	2
9	QE500-07	加压器	Pressurizer	/	1
10	2131700111500	连接杆（直杆）	Connecting rod (Straight rod)	ZIVL01 Φ11/150mm	4
	2131700112000	连接杆（直杆）	Connecting rod (Straight rod)	ZIVL01 Φ11/200mm	4
	2131700112500	连接杆（直杆）	Connecting rod (Straight rod)	ZIVL01 Φ11/250mm	4
	2131700113000	连接杆（直杆）	Connecting rod (Straight rod)	ZIVL01 Φ11/300mm	4
	2131700113500	连接杆（直杆）	Connecting rod (Straight rod)	ZIVL01 Φ11/350mm	4
11	QE500-08	器械盒	Instrument box	/	1

组合式外固定支架ZV系列
COMBINATIONAL EXTERNAL FIXATOR ZV SERIES

针管夹（单槽）
Pin to tube coupling (Single-slot type)

产品编号 Product No.	型号规格 Model&Specification
213530006250	213530006251 ZVP01 Φ6×(Φ2.5 ~ Φ4)
213530008040	213530008041 ZVP01 Φ8×(Φ4 ~ Φ5)
213530008050	213530008051 ZVP01 Φ8×(Φ5 ~ Φ6)




管管夹（单槽）
Tube to tube coupling (Single-slot type)

产品编号 Product No.	型号规格 Model&Specification
213630006060	213630006061 ZVR01 Φ6/Φ6
213630008080	213630008081 ZVR01 Φ8/Φ8




针管夹（双槽）
Pin to tube coupling (Double-slot type)

产品编号 Product No.	型号规格 Model&Specification
213730006250	213730006251 ZVP02 Φ6×(Φ2.5 ~ Φ4)
213730008040	213730008041 ZVP02 Φ8×(Φ4 ~ Φ5)
213730008050	213730008051 ZVP02 Φ8×(Φ5 ~ Φ6)



管管夹（双槽）
Tube to tube coupling (Double-slot type)

产品编号 Product No.	型号规格 Model&Specification
213830006060	213830006061 ZVR02 Φ6/Φ6
213830008080	213830008081 ZVR02 Φ8/Φ8



连接管（直管）
Connecting tube (Straight tube)

产品编号 Product No.	型号规格 Model&Specification
213930006100 ~ 2139300063000	ZVL01 Φ6/10 ~ 300mm (at intervals of 5mm)
213930008100 ~ 2139300084000	ZVL01 Φ8/10 ~ 400mm (at intervals of 5mm)
213930008101 ~ 2139300084001	



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZV系列
COMBINATIONAL EXTERNAL FIXATOR ZV SERIES


连接管（半环管）
Connecting tube (Semi-ring tube)

产品编号 Product No.	型号规格 Model&Specification
214030061850	214030061851 ZVL02 Φ6×185mm
214030081800	214030081801 ZVL02 Φ8×180mm



连接管（弧形管）
Connecting tube (Arc tube)


产品编号 Product No.	型号规格 Model&Specification
214130062850	214130062851 ZVL03 Φ6×285mm
214130082850	214130082851 ZVL03 Φ8×285mm



连接管（加压延长管）
Connecting tube (Compression extension tube)

产品编号:212830083300 212830083301
型号规格:ZVL04 Φ8×330mm

Product No:212830083300 212830083301
Model&Specification:ZVL04 Φ8×330mm



连接管（加压延长管）
Connecting tube (Compression extension tube)

产品编号:214230083300 214230083301
型号规格:ZVL04 Φ8×330mm

Product No:214230083300 214230083301
Model&Specification:ZVL04 Φ8×330mm



指掌骨支架
Mini external fixator

配置 Configuration	数量 Qty
Main part Φ2/Φ2	1
Bone screw Φ2	4



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZV系列
COMBINATIONAL EXTERNAL FIXATOR ZV SERIES

腕关节支架
Wrist joint external fixator

配置 Configuration	数量 Qty
Main part Φ2.5/Φ3.5	1
Bone screw Φ3.5	4




标准半环型支架
Standard semi-ring external fixator

配置 Configuration	数量 Qty
Tube to tube coupling (Single-slot type) ZVR01 Φ8/Φ8	6
Pin to tube coupling (Single-slot type) ZVP01 Φ8×(Φ4 ~ Φ5)	6
Connecting tube (Semi-ring tube) ZVL02 Φ8×180mm	2
Connecting tube (Straight tube) ZVL01 Φ8	1
Connecting tube (Compression extension tube) ZVL04 Φ8×330mm	2
Bone screw Φ4	2
K-wire Φ4	2




骨盆支架
Pelvic external fixator

配置 Configuration	数量 Qty
Tube to tube coupling (Single-slot type) ZVR01 Φ8/Φ8	6
Pin to tube coupling (Single-slot type) ZVP01 Φ8×(Φ4 ~ Φ5)	6
Connecting tube (Arc tube) ZVL03 Φ8×285mm	2
Connecting tube (Compression extension tube) ZVL04 Φ8×330mm	1
Connecting tube (Straight tube) ZVL01 Φ8	2
Bone screw Φ5	6



双管粗隆支架
Double-tube tuberosity external fixator

配置 Configuration	数量 Qty
Pin to tube coupling (Double-slot type) ZVP02 Φ8×(Φ4 ~ Φ5)	5
Connecting tube (Straight tube) ZVL01 Φ8	2
Bone screw Φ5	5



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZV系列
COMBINATIONAL EXTERNAL FIXATOR ZV SERIES

双管直型支架
Double-tube straight external fixator

配置 Configuration	数量 Qty
Pin to tube coupling (Double-slot type) ZVP02 Φ8×(Φ4 ~ Φ5)	4
Connecting tube (Straight tube) ZVL01 Φ8	2
Bone screw Φ5	4



肱骨T型支架
Humeral T type external fixator

配置 Configuration	数量 Qty
Pin to tube coupling (Single-slot type) ZVP01 Φ8×(Φ4 ~ Φ5)	4
Tube to tube coupling (Single-slot type) ZVR01 Φ8/Φ8	1
Connecting tube (Straight tube) ZVL01 Φ8	2
Bone screw Φ5	4



肱骨方框支架
Humeral external fixator

配置 Configuration	数量 Qty
Tube to tube coupling (Single-slot type) ZVR01 Φ8/Φ8	2
Pin to tube coupling (Single-slot type) ZVP01 Φ8×(Φ4 ~ Φ5)	6
Connecting tube (Straight tube) ZVL01 Φ8	3
Bone screw Φ4	6



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZV系列
COMBINATIONAL EXTERNAL FIXATOR ZV SERIES

股骨髁支架
Femoral condylar external fixator

配置 Configuration	数量 Qty
Tube to tube coupling (Single-slot type) ZVR01 Φ8/Φ8	3
Pin to tube coupling (Double-slot type) ZVP02 Φ8×(Φ4 ~ Φ5)	2
Pin to tube coupling (Single-slot type) ZVP01 Φ8×(Φ4 ~ Φ5)	5
Connecting tube (Semi-ring tube) ZVL02 Φ8×180mm	1
Connecting tube (Straight tube) ZVL01 Φ8	2
Connecting tube (Arc tube) ZVL03 Φ8×285mm	1
Bone screw Φ5	7



股骨干支架
Femoral shaft external fixator

配置 Configuration	数量 Qty
Tube to tube coupling (Single-slot type) ZVR01 Φ8/Φ8	3
Pin to tube coupling (Double-slot type) ZVP02 Φ8×(Φ4 ~ Φ5)	4
Pin to tube coupling (Single-slot type) ZVP01 Φ8×(Φ4 ~ Φ5)	3
Connecting tube (Arc tube) ZVL03 Φ8×285mm	2
Connecting tube (Straight tube) ZVL01 Φ8	2
Bone screw Φ5	4
K-wire Φ4	1



踝关节支架
Ankle joint external fixator

配置 Configuration	数量 Qty
Tube to tube coupling (Single-slot type) ZVR01 Φ8/Φ8	6
Pin to tube coupling (Single-slot type) ZVP01 Φ8×(Φ4 ~ Φ5)	8
Connecting tube (Compression extension tube) ZVL04 Φ8×330mm	2
Connecting tube (Straight tube) ZVL01 Φ8	4
Bone screw Φ4	4
K-wire Φ4	2



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架ZV系列
COMBINATIONAL EXTERNAL FIXATOR ZV SERIES

胫骨平台支架
Tibia plateau external fixator

配置 Configuration	数量 Qty
Tube to tube coupling (Single-slot type) ZVR01 Φ8/Φ8	3
Pin to tube coupling (Single-slot type) ZVP01 Φ8×(Φ4 ~ Φ5)	8
Connecting tube (Straight tube) ZVL01 Φ8	1
Connecting tube (Semi-ring tube) ZVL02 Φ8×180mm	1
Connecting tube (Arc tube) ZVL03 Φ8×285mm	1
Bone screw Φ5	6
K-wire Φ4	1



半环T型支架
Humeral semi-ring T type external fixator

配置 Configuration	数量 Qty
Tube to tube coupling (Single-slot type) ZVR01 Φ8/Φ8	3
Pin to tube coupling (Single-slot type) ZVP01 Φ8×(Φ4 ~ Φ5)	5
Connecting tube (Straight tube) ZVL01 Φ8	1
Connecting tube (Semi-ring tube) ZVL02 Φ8×180mm	1
Connecting tube (Arc tube) ZVL03 Φ8×285mm	1
Bone screw Φ4	3
K-wire Φ4	1



膝关节支架
Knee joint external fixator

配置 Configuration	数量 Qty
Tube to tube coupling (Single-slot type) ZVR01 Φ8/Φ8	8
Pin to tube coupling (Single-slot type) ZVP01 Φ8×(Φ4 ~ Φ5)	14
Connecting tube (Semi-ring tube) ZVL02 Φ8×180mm	2
Connecting tube (Arc tube) ZVL03 Φ8×285mm	2
Connecting tube (Compression extension tube) ZVL04 Φ8×330mm	2
Bone screw Φ5	14




备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架H系列
COMBINATIONAL EXTERNAL FIXATOR H SERIES


连接件（六耳环）
Connector (Ring with six ears)

产品编号:214341001550 214341001551
型号规格:H11L01 155
Product No:214341001550 214341001551
Model&Specification:H11L01 155



连接件（三耳环）
Connector (Ring with three ears)

产品编号:214441001550 214441001551
型号规格:H11L02 155
Product No:214441001550 214441001551
Model&Specification:H11L02 155




连接件（半环）
Connector (Semi-ring)

产品编号:214541001550 214541001551
型号规格:H11L03 155
Product No:214541001550 214541001551
Model&Specification:H11L03 155



连接件（足环）
Connector (Foot ring)

产品编号:214641001800 214641001801
型号规格:H11L04 180
Product No:214641001800 214641001801
Model&Specification:H11L04 180



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架H系列
COMBINATIONAL EXTERNAL FIXATOR H SERIES

夹持件（四棱立柱）
Clamp (Post with 4 edges)

产品编号 Product No.	型号规格 Model&Specification
215130000010 215130000011	H11J02 1holes
215130000020 215130000021	H11J02 2holes
215130000030 215130000031	H11J02 3holes
215130000040 215130000041	H11J02 4holes
215130000050 215130000051	H11J02 5holes



夹持件（夹针衬套）
Clamp (Bush)

产品编号 Product No.	型号规格 Model&Specification
215241000040 215241000041	H11J03 Φ4
215241000050 215241000051	H11J03 Φ5
215241000060 215241000061	H11J03 Φ6



加压延长件（标准杆）
Compression extension rod (Standard)

产品编号 Product No.	型号规格 Model&Specification
2153341001400 2153341001401	H11Y01 140
2153341001650 2153341001651	H11Y01 165
2153341001900 2153341001901	H11Y01 190



加压延长件（快装杆）
Compression extension rod (Quick installation)

产品编号 Product No.	型号规格 Model&Specification
2154341001650 2154341001651	H11Y02 165
2154341001900 2154341001901	H11Y02 190
2154341002150 2154341002151	H11Y02 215




备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架H系列
COMBINATIONAL EXTERNAL FIXATOR H SERIES


连接件（全环）
Connector (Full ring)

产品编号 Product No.	型号规格 Model&Specification
214741000800	214741000801 HIL01 Φ80
214741000900	214741000901 HIL01 Φ90
214741001000	214741001001 HIL01 Φ100
214741001100	214741001101 HIL01 Φ110
214741001200	214741001201 HIL01 Φ120
214741001300	214741001301 HIL01 Φ130
214741001400	214741001401 HIL01 Φ140
214741001500	214741001501 HIL01 Φ150
214741001600	214741001601 HIL01 Φ160
214741001700	214741001701 HIL01 Φ170
214741001800	214741001801 HIL01 Φ180
214741001900	214741001901 HIL01 Φ190
214741002000	214741002001 HIL01 Φ200
214741002100	214741002101 HIL01 Φ210
214741002200	214741002201 HIL01 Φ220




连接件（半环）
Connector (Semi-ring)

产品编号 Product No.	型号规格 Model&Specification
214841000800	214841000801 HIL02 Φ80
214841000900	214841000901 HIL02 Φ90
214841001000	214841001001 HIL02 Φ100
214841001100	214841001101 HIL02 Φ110
214841001200	214841001201 HIL02 Φ120
214841001300	214841001301 HIL02 Φ130
214841001400	214841001401 HIL02 Φ140
214841001500	214841001501 HIL02 Φ150
214841001600	214841001601 HIL02 Φ160
214841001700	214841001701 HIL02 Φ170
214841001800	214841001801 HIL02 Φ180
214841001900	214841001901 HIL02 Φ190
214841002000	214841002001 HIL02 Φ200
214841002100	214841002101 HIL02 Φ210
214841002200	214841002201 HIL02 Φ220



连接件（2/3环）
Connector (2/3 ring)

产品编号 Product No.	型号规格 Model&Specification
214941000800	214941000801 HIL03 Φ80
214941000900	214941000901 HIL03 Φ90
214941001000	214941001001 HIL03 Φ100
214941001100	214941001101 HIL03 Φ110
214941001200	214941001201 HIL03 Φ120
214941001300	214941001301 HIL03 Φ130
214941001400	214941001401 HIL03 Φ140
214941001500	214941001501 HIL03 Φ150
214941001600	214941001601 HIL03 Φ160
214941001700	214941001701 HIL03 Φ170
214941001800	214941001801 HIL03 Φ180
214941001900	214941001901 HIL03 Φ190
214941002000	214941002001 HIL03 Φ200
214941002100	214941002101 HIL03 Φ210
214941002200	214941002201 HIL03 Φ220




备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架H系列
COMBINATIONAL EXTERNAL FIXATOR H SERIES

连接件（U型环）
Connector (U-shaped ring)

产品编号 Product No.	型号规格 Model&Specification
215030001200	215030001201 HIL04 Φ120
215030001400	215030001401 HIL04 Φ140
215030001600	215030001601 HIL04 Φ160
215030001800	215030001801 HIL04 Φ180



连接件（6mm丝杆）
Connector (6mm threaded rod)


产品编号 Product No.	型号规格 Model&Specification
215530000100 ~ 215530005000	HIL05 10 ~ 500mm
215530000101 ~ 215530005001	(at intervals of 5mm)



关节件（单关节1型）
Articulated part (Single I)

产品编号:215630000440 215630000441
型号规格:HIG01 44

Product No:215630000440 215630000441
Model&Specification:HIG01 44



关节件（单关节2型）
Articulated part (Single II)

产品编号:215730000420 215730000421
型号规格:HIG02 42

Product No:215730000420 215730000421
Model&Specification:HIG02 42



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

组合式外固定支架H系列
COMBINATIONAL EXTERNAL FIXATOR H SERIES

连接件（L型连接板）
Connector (L connecting plate)

产品编号 Product No.	型号规格 Model&Specification
215830000050	215830000051 HIL07 5holes
215830000060	215830000061 HIL07 6holes
215830000070	215830000071 HIL07 7holes
215830000080	215830000081 HIL07 8holes
215830000090	215830000091 HIL07 9holes
215830000100	215830000101 HIL07 10holes
215830000110	215830000111 HIL07 11holes
215830000120	215830000121 HIL07 12holes



加压延长件（延长套）
Compression extension rod

产品编号 Product No.	型号规格 Model&Specification
2159341000960	2159341000961 HIY01 96
2159341001400	2159341001401 HIY01 140



足踝环形支架
Ankle ring external fixator

主要配置 Key configuration	数量 Qty
Connector (Ring with six ears) HIL01 155	1
Connector (Foot ring) HIL04 180	1
Compression extension rod (Standard) HIY01 140	6



下肢延长环形支架
Lower limbs extension ring external fixator

主要配置 Key configuration	数量 Qty
Connector (Full ring) HIL01 Φ180	4
Compression extension rod HIY01 96	4
Connector (6mm threaded rod) HIL05 150mm	8



备注：产品编号末位数0代表非灭菌，末位数1代表灭菌。 Remarks:The last digit of the product number "0" means non-sterile, and "1" means sterile.

钻头
DRILL BITS

钻头
Drill bit



产品编号 Product No.	直径 Diameter(mm)	长度 Length(mm)	材质 Material
Z10115	Φ1.0	115	不锈钢 Stainless Steel
Z12115	Φ1.2	115	
Z15115	Φ1.5	115	
Z18115	Φ1.8	115	
Z20115	Φ2.0	115	
Z25115	Φ2.5	115	
Z25150	Φ2.5	150	
Z25200	Φ2.5	200	
Z28115	Φ2.8	115	
Z28150	Φ2.8	150	
Z28200	Φ2.8	200	
Z30115	Φ3.0	115	
Z30150	Φ3.0	150	
Z30200	Φ3.0	200	
Z32115	Φ3.2	115	
Z32150	Φ3.2	150	
Z32200	Φ3.2	200	
Z35115	Φ3.5	115	
Z35150	Φ3.5	150	
Z35200	Φ3.5	200	
Z38115	Φ3.8	115	
Z38150	Φ3.8	150	
Z38200	Φ3.8	200	
Z40115	Φ4.0	115	
Z40150	Φ4.0	150	
Z40200	Φ4.0	200	
Z42115	Φ4.2	115	
Z42150	Φ4.2	150	
Z42200	Φ4.2	200	
Z45115	Φ4.5	115	
Z45150	Φ4.5	150	
Z45200	Φ4.5	200	
Z48115	Φ4.8	115	
Z48150	Φ4.8	150	
Z48200	Φ4.8	200	
Z50150	Φ5.0	150	
Z50200	Φ5.0	200	

骨牵引针

BONE SCREWS

骨牵引针

Bone screw



产品编号 Product No.	直径 Diameter(mm)	长度 Length(mm)	材质 Material
G12050T	Φ2.0	50	钛合金 Titanium Alloy
G12060T	Φ2.0	60	
G12070T	Φ2.0	70	
G12050S	Φ2.0	50	不锈钢 Stainless Steel
G12060S	Φ2.0	60	
G12070S	Φ2.0	70	
G12560T	Φ2.5	60	钛合金 Titanium Alloy
G12570T	Φ2.5	70	
G12580T	Φ2.5	80	
G12590T	Φ2.5	90	不锈钢 Stainless Steel
G125100T	Φ2.5	100	
G12560S	Φ2.5	60	
G12570S	Φ2.5	70	不锈钢 Stainless Steel
G12580S	Φ2.5	80	
G12590S	Φ2.5	90	
G125100S	Φ2.5	100	钛合金 Titanium Alloy
G13060T	Φ3.0	60	
G13070T	Φ3.0	70	
G13080T	Φ3.0	80	钛合金 Titanium Alloy
G13090T	Φ3.0	90	
G130100T	Φ3.0	100	
G130110T	Φ3.0	110	不锈钢 Stainless Steel
G130120T	Φ3.0	120	
G13060S	Φ3.0	60	
G13070S	Φ3.0	70	钛合金 Titanium Alloy
G13080S	Φ3.0	80	
G13090S	Φ3.0	90	
G130100S	Φ3.0	100	不锈钢 Stainless Steel
G130110S	Φ3.0	110	
G130120S	Φ3.0	120	
G13560T	Φ3.5	60	钛合金 Titanium Alloy
G13570T	Φ3.5	70	
G13580T	Φ3.5	80	
G13590T	Φ3.5	90	不锈钢 Stainless Steel
G135100T	Φ3.5	100	
G135110T	Φ3.5	110	
G135120T	Φ3.5	120	钛合金 Titanium Alloy
G13560S	Φ3.5	60	
G13570S	Φ3.5	70	
G13580S	Φ3.5	80	不锈钢 Stainless Steel
G13590S	Φ3.5	90	
G135100S	Φ3.5	100	
G135110S	Φ3.5	110	钛合金 Titanium Alloy
G135120S	Φ3.5	120	
G14070T	Φ4.0	70	
G14080T	Φ4.0	80	钛合金 Titanium Alloy
G14090T	Φ4.0	90	
G140100T	Φ4.0	100	
G140110T	Φ4.0	110	

骨牵引针

BONE SCREWS

骨牵引针

Bone screw

产品编号 Product No.	直径 Diameter(mm)	长度 Length(mm)	材质 Material
G140120T	Φ4.0	120	钛合金 Titanium Alloy
G140130T	Φ4.0	130	
G140140T	Φ4.0	140	
G140150T	Φ4.0	150	
G140160T	Φ4.0	160	
G140180T	Φ4.0	180	
G140200T	Φ4.0	200	
G140220T	Φ4.0	220	
G140250T	Φ4.0	250	
G14070S	Φ4.0	70	不锈钢 Stainless Steel
G14080S	Φ4.0	80	
G14090S	Φ4.0	90	
G140100S	Φ4.0	100	
G140110S	Φ4.0	110	
G140120S	Φ4.0	120	
G140130S	Φ4.0	130	
G140140S	Φ4.0	140	
G140150S	Φ4.0	150	
G140160S	Φ4.0	160	钛合金 Titanium Alloy
G140180S	Φ4.0	180	
G140200S	Φ4.0	200	
G140220S	Φ4.0	220	
G140250S	Φ4.0	250	
G15080T	Φ5.0	80	
G15090T	Φ5.0	90	钛合金 Titanium Alloy
G150100T	Φ5.0	100	
G150110T	Φ5.0	110	
G150120T	Φ5.0	120	
G150130T	Φ5.0	130	
G150140T	Φ5.0	140	
G150150T	Φ5.0	150	
G150160T	Φ5.0	160	
G150180T	Φ5.0	180	
G150200T	Φ5.0	200	不锈钢 Stainless Steel
G150220T	Φ5.0	220	
G150250T	Φ5.0	250	
G15080S	Φ5.0	80	
G15090S	Φ5.0	90	
G150100S	Φ5.0	100	
G150110S	Φ5.0	110	
G150120S	Φ5.0	120	
G150130S	Φ5.0	130	
G150140S	Φ5.0	140	
G150150S	Φ5.0	150	钛合金 Titanium Alloy
G150160S	Φ5.0	160	
G150180S	Φ5.0	180	
G150200S	Φ5.0	200	
G150220S	Φ5.0	220	
G150250S	Φ5.0	250	
G16080T	Φ6.0	80	
G16090T	Φ6.0	90	钛合金 Titanium Alloy
G160100T	Φ6.0	100	
G160110T	Φ6.0	110	
G160120T	Φ6.0	120	
G160130T	Φ6.0	130	
G160140T	Φ6.0	140	

骨牵引针

BONE SCREWS

骨牵引针

Bone screw

产品编号 Product No.	直径 Diameter(mm)	长度 Length(mm)	材质 Material
G160150T	Φ6.0	150	钛合金 Titanium Alloy
G160160T	Φ6.0	160	
G160180T	Φ6.0	180	
G160200T	Φ6.0	200	
G160220T	Φ6.0	220	
G160250T	Φ6.0	250	
G16080S	Φ6.0	80	不锈钢 Stainless Steel
G16090S	Φ6.0	90	
G160100S	Φ6.0	100	
G160110S	Φ6.0	110	
G160120S	Φ6.0	120	
G160130S	Φ6.0	130	
G160140S	Φ6.0	140	
G160150S	Φ6.0	150	
G160160S	Φ6.0	160	
G160180S	Φ6.0	180	
G160200S	Φ6.0	200	
G160220S	Φ6.0	220	
G160250S	Φ6.0	250	

自攻自钻骨牵引针

SELF-TAPPING SELF-DRILLING BONE SCREW

自攻自钻骨牵引针

Self-tapping self-drilling bone screw



产品编号 Product No.	直径 Diameter(mm)	长度 Length(mm)	材质 Material
G23060T	Φ3.0	60	钛合金 Titanium Alloy
G23070T	Φ3.0	70	
G23080T	Φ3.0	80	
G23090T	Φ3.0	90	
G230100T	Φ3.0	100	
G230120T	Φ3.0	120	
G23060S	Φ3.0	60	不锈钢 Stainless Steel
G23070S	Φ3.0	70	
G23080S	Φ3.0	80	
G23090S	Φ3.0	90	
G230100S	Φ3.0	100	
G230120S	Φ3.0	120	
G24080T	Φ4.0	80	钛合金 Titanium Alloy
G24090T	Φ4.0	90	
G240100T	Φ4.0	100	
G240120T	Φ4.0	120	
G240140T	Φ4.0	140	
G240150T	Φ4.0	150	
G240160T	Φ4.0	160	不锈钢 Stainless Steel
G240180T	Φ4.0	180	
G24080S	Φ4.0	80	
G24090S	Φ4.0	90	
G240100S	Φ4.0	100	
G240120S	Φ4.0	120	
G240140S	Φ4.0	140	不锈钢 Stainless Steel
G240150S	Φ4.0	150	
G240160S	Φ4.0	160	
G240180S	Φ4.0	180	
G25080T	Φ5.0	80	钛合金 Titanium Alloy
G250100T	Φ5.0	100	
G250120T	Φ5.0	120	
G250130T	Φ5.0	130	
G250150T	Φ5.0	150	
G250180T	Φ5.0	180	
G250200T	Φ5.0	200	不锈钢 Stainless Steel
G25080S	Φ5.0	80	
G250100S	Φ5.0	100	
G250120S	Φ5.0	120	
G250130S	Φ5.0	130	
G250150S	Φ5.0	150	
G250180S	Φ5.0	180	钛合金 Titanium Alloy
G250200S	Φ5.0	200	
G26080T	Φ6.0	80	
G260100T	Φ6.0	100	
G260120T	Φ6.0	120	
G260130T	Φ6.0	130	
G260150T	Φ6.0	150	钛合金 Titanium Alloy
G260180T	Φ6.0	180	
G260200T	Φ6.0	200	

自攻自钻骨牵引针
 SELF-TAPPING SELF-DRILLING BONE SCREW

自攻自钻骨牵引针
 Self-tapping self-drilling bone screw

产品编号 Product No.	直径 Diameter(mm)	长度 Length(mm)	材质 Material
G260250T	Φ6.0	250	钛合金 Titanium Alloy
G26080S	Φ6.0	80	
G260100S	Φ6.0	100	
G260120S	Φ6.0	120	
G260130S	Φ6.0	130	
G260150S	Φ6.0	150	
G260180S	Φ6.0	180	
G260200S	Φ6.0	200	
G260250S	Φ6.0	250	
			不锈钢 Stainless Steel



克氏针
 K-WIRE

克氏针
 K-wire



产品编号 Product No.	直径 Diameter(mm)	长度 Length(mm)	材质 Material
G308100T	Φ0.8	100	钛合金 Titanium Alloy
G308120T	Φ0.8	120	
G308150T	Φ0.8	150	
G308250T	Φ0.8	250	
G308300T	Φ0.8	300	
G308100S	Φ0.8	100	不锈钢 Stainless Steel
G308120S	Φ0.8	120	
G308150S	Φ0.8	150	
G308250S	Φ0.8	250	
G308300S	Φ0.8	300	
G310100T	Φ1.0	100	钛合金 Titanium Alloy
G310120T	Φ1.0	120	
G310150T	Φ1.0	150	
G310250T	Φ1.0	250	
G310300T	Φ1.0	300	
G310100S	Φ1.0	100	不锈钢 Stainless Steel
G310120S	Φ1.0	120	
G310150S	Φ1.0	150	
G310250S	Φ1.0	250	
G310300S	Φ1.0	300	
G312100T	Φ1.2	100	钛合金 Titanium Alloy
G312120T	Φ1.2	120	
G312150T	Φ1.2	150	
G312250T	Φ1.2	250	
G312300T	Φ1.2	300	
G312100S	Φ1.2	100	不锈钢 Stainless Steel
G312120S	Φ1.2	120	
G312150S	Φ1.2	150	
G312250S	Φ1.2	250	
G312300S	Φ1.2	300	
G315100T	Φ1.5	100	钛合金 Titanium Alloy
G315120T	Φ1.5	120	
G315150T	Φ1.5	150	
G315250T	Φ1.5	250	
G315300T	Φ1.5	300	
G315100S	Φ1.5	100	不锈钢 Stainless Steel
G315120S	Φ1.5	120	
G315150S	Φ1.5	150	
G315250S	Φ1.5	250	
G315300S	Φ1.5	300	
G316250T	Φ1.6	250	钛合金 Titanium Alloy
G316300T	Φ1.6	300	不锈钢 Stainless Steel
G316250S	Φ1.6	250	
G316300S	Φ1.6	300	钛合金 Titanium Alloy
G318250T	Φ1.8	250	
G318300T	Φ1.8	300	不锈钢 Stainless Steel
G318250S	Φ1.8	250	
G318300S	Φ1.8	300	钛合金 Titanium Alloy
G320250T	Φ2.0	250	

克氏针

K-wires

克氏针

K-wire

产品编号 Product No.	直径 Diameter(mm)	长度 Length(mm)	材质 Material
G320300T	Φ2.0	300	钛合金 Titanium Alloy
G320250S	Φ2.0	250	不锈钢 Stainless Steel
G320300S	Φ2.0	300	
G325250T	Φ2.5	250	钛合金 Titanium Alloy
G325300T	Φ2.5	300	不锈钢 Stainless Steel
G325250S	Φ2.5	250	
G325300S	Φ2.5	300	
G330250T	Φ3.0	250	钛合金 Titanium Alloy
G330300T	Φ3.0	300	不锈钢 Stainless Steel
G330250S	Φ3.0	250	
G330300S	Φ3.0	300	
G335250T	Φ3.5	250	钛合金 Titanium Alloy
G335300T	Φ3.5	300	不锈钢 Stainless Steel
G335250S	Φ3.5	250	
G335300S	Φ3.5	300	
G340250T	Φ4.0	250	钛合金 Titanium Alloy
G340300T	Φ4.0	300	不锈钢 Stainless Steel
G340250S	Φ4.0	250	
G340300S	Φ4.0	300	
G345250T	Φ4.5	250	钛合金 Titanium Alloy
G345300T	Φ4.5	300	不锈钢 Stainless Steel
G345250S	Φ4.5	250	
G345300S	Φ4.5	300	